

- SENSITIVE SKIN LIFTING
- EYE AREA - EYELIDS
- LIFTING
- SKIN TURGOR
- ARM TREATMENT
- THIGH TIGHTENING
- STOMACH TIGHTENING



## FACE & BODY TREATMENT

# RF DYNAMIC FACE & BODY

RADIO WAVE FACE AND BODY LIFTING

Code: 5058

RF DYNAMIC FACE & BODY is a technologically advanced device that offers extreme ease of use and great versatility. It ensures instant rejuvenation and tightening of facial skin and visibly reduced cellulite on the body.

Through a sophisticated electronic system, it transfers heat deeper into the fatty and subcutaneous tissue without damaging the epidermis, giving the skin a younger and tighter look.

The frequency of treatment enables better blood circulation, restoring the tension in the facial muscles. Facial treatment is also recommended for elderly clients with visible skin imperfections. Radio waves by means of their thermal effect cause the synthesis of new collagen and elastin, while the corrected spirals of already present collagen contract, which leads to instant skin tightening and the lifting effect is seen immediately after the first treatment.

RF Dynamic treatment is intended for face and body as it improves wrinkles and fine lines around the eyes, drooping eyelids, eye bags, fine forehead and face wrinkles, sagging cheeks, double chin, neck, décolleté, chest, arms, abdomen, thighs and knees.

## FRACTIONAL METHOD OF REJUVENATION

The RF DYNAMIC operating system is fractionated energy - divided into the smallest parts, where, in addition to the renewal of the surface part, the heat is also applied in deeper layers of tissues.

The effect that is achieved in this way is a complete renewal of the outer surface of the skin, but also a strong stimulation of the production of new collagen in the deeper layers of the skin. As collagen is the basic ingredient that makes the skin firm and taut, we achieve rejuvenation with a maintained natural appearance after this procedure, using the natural potentials of our body.

## FACE

Radio wave lifting is performed for atonic skin, premature aging, photoaging, reduced skin turgor, before and post-surgical treatments. Smaller wrinkles on the forehead and around the eyes disappear and the deep ones become visibly shallower and less prominent, the nasolabial folds are reduced and the skin tightens.

## BODY

Radio wave body lifting can be done after the cavitation treatment, for tightening the chest area, gluteus, sagging skin (abdomen, inner thighs and arms), in order to reduce stretch marks and as an anti-cellulite treatment.

## ADDITIONAL EQUIPMENT



Code: 8090



Code: 8091

## Technical characteristics

Product name: RF DYNAMIC FACE & BODY

Code: 5058

Energy source: radiofrequency

Device power: max 50 W

Large lever output power: 160 V

Small lever output power: 110 V

RF frequency: 1 MHz

Working temperature: from -10 do 40°C

Voltage: 220-240 V, 50-60 Hz

Dimensions (LxWxH) 320x400x200 mm

Weight: 5 kg