


STANDARD SET

- **ULTRAPEEL** - facial peeling with an ultrasonic spatula
- **DIAMOND PEEL** - diamond microdermabrasion
- **VACUUM EXPEDITE** - improving peripheral circulation and drainage
- **MESOTHERAPY** - removal and reduction of changes
- **RF LIFT** - radio wave lifting
- **ULTRA FACE** - increasing the rate of metabolic reactions
- **LASER REGENERATION** - regeneration & renewal of the dermis structure
- **HIGH FREQUENCY** - draining effect

ADDITIONAL FUNCTIONS

(NOT INCLUDED IN THE STANDARD KIT):

- **OZONATION - VAPOZON** - the beneficial effect of water vapor ensures an effective antibacterial result
- **COSMETIC MAGNIFIER** - to enlarge the work area



Code 9024

FACE TREATMENT

FACE CREATOR

MULTIFUNCTIONAL DEVICES

IT IS A MULTIFUNCTIONAL DEVICE THAT COMBINES 8 METHODS FOR TREATMENTS OF REGENERATION, PEELING AND REJUVENATION. IT IS DIVIDED INTO SEGMENTS: ULTRASONIC PEELING - EXFOLIATION, DIAMOND MICRODERMABRASION, VACUUM, MESOPORATION, RADIO WAVE LIFTING, ULTRASONIC TREATMENT - SONOPHORESIS, LASER SKIN REGENERATION (BIOLASER), DISINFECTION AND IMPROVEMENT OF PERIPHERAL CIRCULATION - HIGH FREQUENCY.

ULTRAPEEL- 27kHz - ultrasonic peeling and facial cleansing

This treatment provides deep cleaning and insertion of active ingredients. It combines the effects of ultrasonic vibration using a spatula for even facial peeling. The use of various frequencies allows for regeneration, toning, firming and they ensure favourable results on all skin types by means of 4 specific programs.

DIAMOND PEEL- diamond dermabrasion

This treatment provides even skin peeling as well as stimulating cell renewal in the treated tissue. It increases skin permeability and provides better absorption of cosmetic preparations.

VACUUM EXPEDITE- vacuum treatment

This treatment is pleasant and as effective as a manual massage; it allows for a more suitable facial treatment and it is also less fatiguing for the beautician. It uses the principle of vacuum massage with mild aspiration to prevent damage to the dermis and the capillary system.

MESOTHERAPY- electroporation

Through electrical impulses, electroporation contributes to the opening of the channel in the dermis in order to transfer deeply the substance being applied, in accordance with the targeted irregularities.

RF LIFT- radio wave lifting

This treatment uses the process of diathermy, which allows vasodilation and increases the flow of oxygen to the tissues, accelerates cellular metabolism and stimulates the production of collagen, elastin and hyaluronic acid, making the skin elastic and firm.

ULTRA FACE- ultrasonic facial treatment

This ULTRASOUND type is used at a frequency of 1MHz in facial treatments for the toning of sagging and tired-looking skin, for redness (couperosis) treatment, in order to regulate sebum secretion and in the treatment of blemishes caused by increased pigmentation.



HIGH FREQUENCY- high frequency current

The application of high-frequency current stimulates the surface circulation of tissues. The formation of ozone, which has an antibacterial and drying effect on the skin, also contributes to this. The main application is related to acne, seborrheic skin and it also ensures stimulating effects through facial preparations.

LASER REGENERATION - biostimulation laser

BIOLASER 100mW is a device that stimulates collagen synthesis, thus accelerating wound healing, relieving pain (analgesic effect), reducing swelling and scarring; it provides both a biological effect and a beneficial effect on the immune system.

VAPOZON (Option - not included in the basic set of the device)

Vapozon is a device for dispersing water vapor which hydrates the skin and accelerates its respiration and metabolism in cells. During the treatment, ozone is mixed with water vapor and ionized. Ionized water vapor has a fine structure which allows optimal hydration of the skin.



Technical characteristics

Product name: FACE CREATOR
 Code: 9024
 Voltage: 220-240 V AC, 50/60 Hz
 Housing: Steel
 Power: 150 W
 Main fuse: T-2 A
 IC Classification IEC60601/1: Class I, Type BF
 Screen dimensions: 10,2" color touch screen
 Number of handles: 8
 Number of programs: 15
 Dimensions (LxWxH) 620x370x1110 mm
 Weight: 30 kg