



FACE TREATMENT

MEZO RF LIFT TOUCH UG

MEZOPORATION AND RADIO WAVE LIFTING

Code: 9870

RADIOFREQUENCY

This treatment uses the process of diathermy, which allows vasodilation and increases the flow of oxygen to the tissues, accelerates cellular metabolism and stimulates the production of collagen, elastin and hyaluronic acid, making the skin elastic and firm. Depending on the problem, the most beneficial results are achieved by combining complementary techniques (treatments), as multiple techniques increase optimal results. This primarily refers to the combination of mesoporation and radiowave lifting, where not only are the signs of aging and skin deterioration removed, but also the aging process is slowed down, which is the ultimate goal when we talk about non-surgical methods of rejuvenation.

MESOPORATION

This is a treatment in which the technique of introducing active substances into the deeper layers of the skin is applied by means of a system of adapted current impulses which act quickly, efficiently and in a targeted way to repair the skin structure.

SCREEN DISPLAY



Technical characteristics

Product name: MEZO RF LIFT TOUCH UG
 Code: 9870
 Voltage: 220-240 V AC; 50/60 Hz
 Device power: 50 W
 Energy source: Radiofrequency and DC pulse voltage
 Radiofrequency: 1 MHz
 Number of handles: 2
 Number of programs: 5
 Pulse shape: square, 500 ms
 Main fuse: T-1 A
 IC Classification IEC60601/1: CLASS I, TYPE BF
 housing: polystirol
 Dimensions (LxWxH): 320x400x200 mm
 Weight: 8 kg