

- FACE LIFTING
- EYE ZONE- EYELIDS
- SENSITIVE SKIN LIFTING
- SKIN TURGOR



FACE TREATMENT

RF DYNAMIC FACE

RADIO WAVE FACELIFT

Code: 5062

PRINCIPLE OF OPERATION

RF DYNAMIC is a device that uses the radio frequency on the face area in order to provide effective treatments of the irregularities present at different depths of the skin.

Through a sophisticated electronic system, RF DYNAMIC FACE transfers heat to the deeper layers of the skin without damaging the epidermis, providing a younger and more toned look to the skin.

This method is used for aesthetic problems such as facial wrinkles, drooping eyelids, bags under the eyes, double chin and neck wrinkles. Sagging skin of the chest area is also treated.

The heat produced causes an increase in the amount of oxygen that is transferred to the tissues, accelerates cell metabolism and stimulates the production of collagen, elastin and hyaluronic acid, providing a lifting effect and progressive skin tightening.

EFFECTS

With the RF DYNAMIC FACE operating system, it is possible to perform facial treatments that provide visible results in relation to the wrinkles and the entire surface of the face where sagging skin is apparent.

The frequency of the heat wave release ensures improved blood circulation, restoring facial energy.

The treatment is pleasant and warming and it is recommended for people with initial or visible signs of a sagging face, wrinkles and decreased facial muscle tone.

ADVANTAGES

- > Energy 50w / cm²
- > Multipolar probe providing uniform heat in deeper tissue layers
- > Defined program procedures for the face area

ADDITIONAL EQUIPMENT



Code: 8091

Technical characteristics:

Product name: RF DYNAMIC FACE
 Code: 5062
 Energy source: radio frequency
 Output power: max 50 W
 Small lever output power: 110 V
 RF frequency: 1 MHz
 Operating temperature: -10 do 40°C
 Voltage: 220-240V, 50-60 Hz
 Dimensions (LxWxH) 320x400x200 mm
 Weight: 4 kg